



## PERFECT ANEROID CLINIC 48

### PERFECTLY ACCURATE AND RELIABLE

Clinically validated, 100% inspected and with an accuracy of 1% over the entire measuring range & precision defined by ERKA. Hand-polished with a bright sheen and chrome-plated, the brass housing provides maximum surface hygiene. The ergonomically optimized design as well as the flat design ensure easy handling. The result: perfection in blood pressure monitoring! Based on this claim, we have developed the clinic model with manometer protective cap for use in the hectic everyday clinical environment.



| Perfect Aneroid Clinic 48           | Size   | Reference  |
|-------------------------------------|--------|------------|
| with GreenCuff Superb Rapid , green | Size 4 | 215.20482. |
| with GreenCuff Superb Rapid , green | Size 5 | 215.24482. |
| with GreenCuff Superb Rapid , green | Size 3 | 215.26482. |
| with GreenCuff Superb Rapid , green | Size 2 | 215.27482. |
| with GreenCuff Superb Rapid , green | Size 1 | 215.28482. |
| with GreenCuff Superb Rapid , green | Size 6 | 215.29482. |

| Accessories                          | Size | Reference |
|--------------------------------------|------|-----------|
| Cuff: Green Cuff Superb Hook, green  | 5    | 025.24282 |
| Cuff: Green Cuff Superb Hook, green  | 4    | 025.20282 |
| Cuff: Green Cuff Superb Hook, green  | 3    | 025.26282 |
| Cuff: Green Cuff Superb Rapid, green | 6    | 025.29482 |
| Cuff: Green Cuff Superb Rapid, green | 5    | 025.24482 |
| Cuff: Green Cuff Superb Rapid, green | 4    | 025.20482 |
| Cuff: Green Cuff Superb Rapid, green | 3    | 025.26482 |
| Cuff: Green Cuff Superb Rapid, green | 2    | 025.27482 |
| Cuff: Green Cuff Superb Rapid, green | 1    | 025.28482 |

Spare Parts

Pump Ball green

000.20200

Legal Notes and Data Protection

ERKA.  
Kallmeyer Medizintechnik GmbH & Co. KG  
Im Farchet 15  
83646 Bad Tölz

Geschäftsführer: Björn Kallmeyer

HRA 75980, Amtsgericht München  
Umsatzsteuer-ID: DE 812 833 893

Kontakt:  
Telefon: 0 80 41 / 80 09-0  
Telefax: 0 80 41 / 80 09-39  
eMail: [made@erka.org](mailto:made@erka.org)  
<http://www.erka.org>

Verantwortlich:  
Björn Kallmeyer

Signatur / Stempel

PERFECT ANEROID CLINIC 48

