



SMART RAPID

ALWAYS RIGHT

Correct blood pressure measurement is the cornerstone of the diagnosis, management and treatment of arterial hypertension. Selection of the correct cuff size is of key significance. So-called undercuffing or overcuffing, i.e. the selection of a cuff which is too small or too large relative to the upper arm circumference can lead to a distortion of the measured values by up to 30 mmHg. Accordingly, ERKA. offers a wide range of different cuff sizes - from infant to thigh - the right size for each user. And what's more: cuff change in the blink of an eye - the ERKA. quick closure makes it possible. This is our SMART Concept. Biocompatible, water-resistant and entirely dermatologically safe: our GREEN CUFF Smart cuffs can be superbly cleaned and easily disinfected. This means maximum hygiene for patient and staff.

Smart RapidSize	Reference
Size 4	024.20493.
Size 5	024.24493.
Size 3	024.26493.
Size 2	024.27493.
Size 1	024.28493.
Size 6	024.29493.

Accessories

Size

Reference

Spare Parts

Legal Notes and Data Protection

ERKA.
Kallmeyer Medizintechnik GmbH & Co. KG
Im Farchet 15
83646 Bad Tölz

Geschäftsführer: Björn Kallmeyer

HRA 75980, Amtsgericht München
Umsatzsteuer-ID: DE 812 833 893

Kontakt:
Telefon: 0 80 41 / 80 09-0
Telefax: 0 80 41 / 80 09-39
eMail: made@erka.org
<http://www.erka.org>

Verantwortlich:
Björn Kallmeyer

Signatur / Stempel

SMART RAPID

