



## SMART RAPID

### ALWAYS RIGHT

Correct blood pressure measurement is the cornerstone of the diagnosis, management and treatment of arterial hypertension. Selection of the correct cuff size is of key significance. So-called undercuffing or overcuffing, i.e. the selection of a cuff which is too small or too large relative to the upper arm circumference can lead to a distortion of the measured values by up to 30 mmHg. Accordingly, ERKA. offers a wide range of different cuff sizes - from infant to thigh - the right size for each user. And what's more: cuff change in the blink of an eye - the ERKA. quick closure makes it possible. This is our SMART Concept. Biocompatible, water-resistant and entirely dermatologically safe: our GREEN CUFF Smart cuffs can be superbly cleaned and easily disinfected. This means maximum hygiene for patient and staff.



Smart RapidSize	Reference
Size 4	024.20493.
Size 5	024.24493.
Size 3	024.26493.
Size 2	024.27493.
Size 1	024.28493.
Size 6	024.29493.

Accessories

Size

Reference

Spare Parts

Legal Notes and Data Protection

ERKA.  
Kallmeyer Medizintechnik GmbH & Co. KG  
Im Farchet 15  
83646 Bad Tölz

Geschäftsführer: Björn Kallmeyer

HRA 75980, Amtsgericht München  
Umsatzsteuer-ID: DE 812 833 893

Kontakt:  
Telefon: 0 80 41 / 80 09-0  
Telefax: 0 80 41 / 80 09-39  
eMail: [made@erka.org](mailto:made@erka.org)  
<http://www.erka.org>

Verantwortlich:  
Björn Kallmeyer

Signatur / Stempel

SMART RAPID

